

SALERA / SALERA PRESELECT							
article	material joint upper part	material joint upper part	pelvic band (in mm)	Joint carrier made of aluminum (thickness x width x length*, in mm)	body weight (in kg)	height (in cm)	Item number
SALERA extra large	aluminium	aluminium	10 x 50 x 600	8 x 22 x 480	up to 90	up to 185	30.05.500.001
SALERA large	stole	titanium	10 x 50 x 600	6 x 19 x 320	up to 55	up to 160	30.05.400.001
SALERA preselect large	stole	titanium	10 x 50 x 600	6 x 19 x 320	up to 55	up to 160	30.05.400.002
SALERA medium	aluminium	titanium	8 x 50 x 500	6 x 19 x 320	up to 40	up to 145	30.05.300.001
SALERA preselect medium	aluminium	titanium	8 x 50 x 500	6 x 19 x 320	up to 40	up to 145	30.05.300.002
SALERA small	aluminium	titanium	6 x 50 x 400	6 x 16 x 270	up to 20	up to 120	30.05.200.001
SALERA preselect small	aluminium	titanium	6 x 50 x 400	6 x 16 x 270	up to 20	up to 120	30.05.200.002

* Length measured from the center of the joint screw (joint center) up to the end of the joint carrier.

GOTTINGER HANDELSHAUS OHG

Ilching Weg 1 ▪ 85604 Zorneding ▪ phone: 08106/3663-36 ▪ Fax: 08106/3663-31 ▪ Email: handelshaus@gottinger.de